

# Biotechnology Explorer™


CAPTIVATING SCIENCE EDUCATION

A scientist wearing a white lab coat and blue gloves is shown in a close-up, holding a test tube. The test tube contains a yellow and green liquid. Overlaid on the scientist's face is a semi-transparent image of a DNA gel electrophoresis pattern, with the sequence 't g a g t t g c c t t e g g a t t' visible. The background is dark, and the overall image has a scientific and educational theme.

**BIO-RAD**


**Professional Development**

## Crime Scene Investigator PCR Basics™



# Crime Scene Investigator PCR Basics™ Kit

## Instructors



### **Stan Hitomi**

Coordinator – Math & Science  
San Ramon Valley Unified School District  
Danville, CA

### **Kirk Brown**

Lead Instructor, Edward Teller Education Center  
Science Chair, Tracy High School  
and Delta College, Tracy, CA

### **Sherri Andrews, Ph.D.**


Curriculum and Training Specialist  
Bio-Rad Laboratories

### **Essy Levy, M.Sc.**

Curriculum and Training Specialist  
Bio-Rad Laboratories

## **Why Teach Crime Scene Investigator PCR Basics™ Kit ?**

- **Exciting real-world connections**
- **Tangible results**
- **Statistical Analysis**
- **Standards-based**



## Target Audience

- **The Crime Scene Investigator PCR Basics™ Kit is intended to be an introduction to the polymerase chain reaction (PCR)**
- **Students will have a much better appreciation of the kit if they have some understanding of DNA structure and function**

## **Crime Scene Investigator PCR Basics™ Kit Advantages**

- **Standards Based**
- **Perform real-world DNA profiling**
- **Use PCR to amplify multiple DNA samples**
- **Use electrophoresis to visualize results**
- **Complete in two 45 minute lab sessions**
- **Sufficient materials for 8 student workstations**

## Workshop Time Line

- **Introduction to DNA profiling**
- **Set up PCR reactions**
- **Electrophorese PCR products**
- **Analysis and interpretation of results**



## What is DNA profiling?

**DNA profiling is the use of molecular genetic methods to determine the exact genotype of a DNA sample in a way that can basically distinguish one human being from another**

**The unique genotype of each sample is called a **DNA profile**.**


## How do crime scene investigators create a DNA profile?

### 1. Evidence is collected at the crime scene:



## How do crime scene investigators create a DNA profile?


**2. DNA is extracted from sources at the crime scene and from victim and suspects**



## How do crime scene investigators create a DNA profile?

### 3. DNA samples are processed

Sample Obtained from  
Crime Scene or  
Paternity Investigation



**Since humans are 99.9% identical where do crime scene investigators look for differences in DNA profiles?**

**4. Crime Scene Investigators search in areas of the genome that are unique from individual to individual and are “anonymous” (control no known trait or function)  
The areas examined are Short Tandem Repeats or STR’s**



**STR region**

## Example of an STR: TH01

**The TH01 locus contains repeats of TCAT.**

**CCC TCAT TCAT TCAT TCAT TCAT TCAT AAA**

**This example has 6 TCAT repeats.**

**There are more than 20 known TH01 alleles.**

**Each individual inherits 1 allele from each parent.**

## Determining genotypes for individuals using STRs

**Ms. Smith's TH01 locus for her two chromosomes is given below.**

**What is her genotype?**

**MOM'S CHROMOSOME**

**CCC TCAT TCAT TCAT TCAT TCAT TCAT AAA**

**DAD'S CHROMOSOME**


**CCC TCAT TCAT TCAT TCAT TCAT TCAT TCAT  
TCAT TCAT TCAT TCAT TCAT TCAT TCAT AAA**

To determine the genotype (DNA profile) Crime Scene Investigators make billions of copies of the target sequence using PCR

## Target DNA



STR region



Starting DNA  
Template

**PCR**



## What's the point of PCR?

- **PCR, or the polymerase chain reaction, makes copies of a specific piece of DNA**
- **PCR allows you to look at one specific piece of DNA by making copies of \*only\* that piece of DNA**
- **PCR is like looking for a needle in a haystack, and then making a haystack out of the needle**


**DNA profiling is used to determine which suspect can not be excluded from suspicion.**



## How are suspects included or excluded from an investigation?

- **Suspects are included in an investigation if their DNA profile matches with genotypes found at the crime scene**
- **Suspects can be excluded if their DNA profile does not match genotypes found at the crime scene**

# Crime Scene Investigator PCR Basics™ Procedures Overview




# Laboratory Quick Guide

## Quick Guide

### Lesson 1: Setting up the PCR Reactions


1. Label 5 PCR tubes CS, A, B, C, or D, and include your group name or initials as well. Place the PCR tube in a capless tube in the foam float on ice.



2. Using the chart below as a guide, transfer 20  $\mu$ l of the appropriate template DNA into the correctly labeled tube. Important: use a fresh aerosol barrier pipet tip for each DNA sample.

Tube label	DNA	Master mix + primers
CS + your initials	20 $\mu$ l Crime Scene DNA	20 $\mu$ l MMP (blue)
A + your initials	20 $\mu$ l Suspect A DNA	20 $\mu$ l MMP (blue)
B + your initials	20 $\mu$ l Suspect B DNA	20 $\mu$ l MMP (blue)
C + your initials	20 $\mu$ l Suspect C DNA	20 $\mu$ l MMP (blue)
D + your initials	20 $\mu$ l Suspect D DNA	20 $\mu$ l MMP (blue)

3. Transfer 20  $\mu$ l of the blue MMP into each of the 5 PCR tubes containing template DNA. Pipet up and down to mix. Cap each tube after adding blue MM. Important: use a fresh aerosol barrier pipet tip each time. Immediately cap each tube after adding MMP.




4. Place your capped PCR tubes in their adaptors on ice.

5. When instructed to do so, place your tubes in the thermal cycler. Your instructor will program it for PCR.



## Set up PCR reactions




1. Find the PCR tubes at your station. Label them 'CS' for **C**rime **S**cene DNA, 'A' for Suspect **A** DNA, 'B' for Suspect **B** DNA, 'C' for Suspect **C** DNA, and 'D' for Suspect **D** DNA.
2. Keeping the tubes on ice, add 20  $\mu$ l of Master Mix + blue primers to each tube.
3. Keeping the tubes on ice, add 20  $\mu$ l of each DNA to the appropriately labeled tube.
4. **USE A FRESH TIP EACH TIME!**
5. Mix and put in thermal cycler
6. Cycle ~3 hours

# The PCR Reaction

## What do you need?

- **Template (containing the STR you want to amplify for the study)**
- **Sequence-specific primers flanking the target sequence**



- **Nucleotides (dATP, dCTP, dGTP, dTTP)**
- **Magnesium chloride (enzyme cofactor)**
- **Buffer, containing salt**
- ***Taq* polymerase**

**What is  
happening in  
the PCR tube  
while in the  
thermocycler?**

## **PCR Animation**

[http://www.bio-rad.com/flash/07-0335/07-0335\\_PCR.html](http://www.bio-rad.com/flash/07-0335/07-0335_PCR.html)



## The PCR Reaction

### How does it work?

Heat (**94°C**) to denature DNA strands

Cool (**52°C**) to anneal primers to template


Warm (**72°C**) to activate *Taq* polymerase, which extends primers and replicates DNA

Repeat **35** cycles



**Heat causes DNA  
strands to separate**

**Denaturation of  
DNA at 94°C**




**Primers bind to the template sequence**




**Taq polymerase recognizes 3' end of primer + template strand**

**Primers anneal at 52°C**




**Taq extends at 72°C**




**STR DNA is replicated**

**Repeat denaturing, annealing, and extending 35 cycles**

**The exact-length target product is made in the third cycle**



To visualize  
PCR products  
Crime Scene  
investigators  
use gel  
electrophoresis



## Electrophorese PCR products

- 1. Add 10 ul of Orange G Loading Dye to each PCR tube and mix**
- 2. Set up gel and electrophoresis equipment**
- 3. Load 20 ul of CSI allele ladder to Lane 1**
- 4. Load 20 ul of your PCR reactions in lanes 2 to 6**
- 5. Electrophorese samples**
- 6. Stain gel with Fast Blast DNA Stain**
- 7. Analyze results**

## Using the digital micropipet


Add 10ul of loading dye to each microtube



## Agarose Electrophoresis

**Place gel in gel  
box**


**Pour buffer in  
box until gel  
wells are  
covered.**





**Place 20ul of samples into appropriate wells**


**Set up electrophoresis chamber by putting top in place and connecting it to the power supply**



## Agarose Electrophoresis Running

**Agarose gel sieves**  
DNA fragments  
according to size  
– Small fragments  
move farther than  
large fragments

Use a 3% gel to  
separate small  
fragment sizes



## Milestones in Forensic DNA analysis

- |              |  |
|--------------|--|
| <b>1985</b>  | Alec Jeffries develops RFLP  |
| <b>1990</b>  | PCR analysis using single locus STR begins   |
| <b>1992</b>  | FBI initiates STR work   |
| <b>1994</b>  | DNA Identification Act: provides funding for national DNA database   |
| <b>1995</b>  | OJ Simpson trial focuses public attention on DNA evidence  |
| <b>1998</b>  | FBI starts CODIS database; Swissair disaster – all remains identified using STR DNA profiling                |
| <b>2001</b>  | World Trade Center disaster in NYC – many remains identified using a combination of DNA profiling approaches |
| <b>2004</b>  | Indian Ocean tsunami; Interpol and other world agencies use DNA profiling to identify victims                |
| <b>Today</b> | Trace your Genetic Genealogy; commercially available packages can trace paternal/maternal ancestry           |

# DNA Testing Today

**GeneTree.com &  
Ancestry.com  
provide DNA tests  
from \$99-\$200 to  
trace genealogy**

TUESDAY, OCTOBER 23, 2007

☆☆☆☆

San Francisco Chronicle B5

## Sites offer DNA tests, databases to assist with genealogy research

By Anick Jesdanun  
ASSOCIATED PRESS

NEW YORK — Two services launching just a week apart tap a growing interest in DNA testing to help people find their ancestors and learn more about their lives.

GeneTree, which opens today, and Ancestry.com, which started its DNA Ancestry service last Tuesday, both sell DNA kits for less than \$200. Users can build online family trees and contact others with DNA matches to compare family histories.

Genealogy research has become popular in recent years as online services improve access to vast databases of immigration, military and other records from around the world. According to the Pew Internet and American Life Project, a quarter of Internet users have researched their ancestors online.

Lately, many of them have been turning to DNA testing to uncover additional clues, said Dick Eastman, who writes a newsletter about online genealogy. Although DNA won't provide all the answers, such as names and precise dates, he said, it could open leads.

"Anybody who's got a mystery is going to do this sooner or later, and that's a pretty high percentage of us," Eastman said.

That's particularly true of black Americans, many of whom have trouble tracing roots beyond the slavery era, Eastman said. Eastern Europeans, Jews and certain other groups also find records fragmented, he said.

GeneTree and Ancestry join services from Family Tree DNA and others.

James Lee Sorenson, GeneTree's chief executive, said he believes his site stands out for its exclusive access to records from the nonprofit Sorenson Molecular Genealogy Foundation. The group has collected DNA samples from 100,000 individuals worldwide and conducted ancestry research on them to produce a larger database of 6 million people.

Ancestry, based in Provo, Utah, is building its DNA database largely from scratch; company officials say they are on track to capture the genetic profiles of 50,000 people within six months.

Both GeneTree and Ancestry use DNA test kits from Sorenson Genomics.

GeneTree sends mouthwash

that users swirl in their mouths, spit into a container and mail back for \$99 or \$149, depending on how much DNA the user wants analyzed. With Ancestry, users return a cheek swab. Ancestry offers a greater variety of tests; the one comparable to GeneTree's high-end offering costs \$179.

Besides finding matches, DNA patterns can help assess the likely origins of an individual's ancestors thousands of years ago, allowing the user to then visually trace migration backward to the first humans, widely believed to hail from Africa.

Both sites are incorporating elements of social networking, akin to those at Facebook and MySpace.

With GeneTree, each family member in the online tree — whether a user of GeneTree or not — gets a personal profile page. Users may add photos, video and other documents to their own pages or those of relatives, using free tools from sister company Sorenson Media.

Ancestry already has some tools for adding photos and other files and plans additional features, such as letting users with the same last name DNA results.



# Genealogical Analysis Uses:

DNA sequencing (mtDNA)

Single Nucleotide Polymorphism testing (SNPs)

Short Tandem Repeat testing (STRs)

## PATERNAL LINEAGE TEST RESULTS

[DNA Home](#)

**John Harold Doe**

[View Maternal Results](#)

After reviewing your results, you can begin to compare your profile against other DNA Ancestry participants. Click the "Find Paternal Matches" button to begin.

[Find Paternal Matches](#)

Paternal Lineage Test (Y-Chromosome 46)  
[Understanding your results](#)

[PDF](#) Download or print a copy of your results.

DYS19a	DYS19b	DYS385a	DYS385b	DYS208	DYS389I	DYS389II	DYS290	DYS291	DYS292	DYS293	DYS426
-	-	-	-	7	9	21	19	6	7	10	6.2
DYS427	DYS438	DYS439	DYS441	DYS442	DYS444	DYS445	DYS446	DYS447	DYS448	DYS449	DYS452
12	7	8	7	8	9	6	8	14	17	23	24
DYS454	DYS455	DYS456	DYS458	DYS459a	DYS459b	DYS460	DYS461	DYS462	DYS463	DYS464a	DYS464b
6	7	11	12	-	-	7	7	8	16	-	-
DYS464c	DYS464d	DYS464e	DYS464f	00AAT1807	YCAIIa	YCAIIb	Y-GATA-A1C	DYS635	Y-GATA-H4		
-	-	-	-	6	-	-	9	15	8		

### PATERNAL ANCIENT ANCESTRY



## R1b


Paternal Ancient Ancestry Haplogroup R1b first arrived in Europe from West Asia during the Upper Paleolithic period (35,000–40,000 years ago) at the beginning of the **Aurignacian** culture. This culture is one of the first within Europe to leave cave-art and their stone tools were more refined than previous periods. The **Magdalenian** culture is also considered by some to have existed at this time.

As the last ice-age began, it became necessary to move down to below the tree-line to hunt game. At its peak, the ice sheet within Europe extended down as far as southern Ireland, the middle of England and across northern Germany. Scandinavia was entirely covered. The sea-ice

**Crime scene investigators use techniques that are fast, cost effective, and have a high Power of Discrimination**

**The Power of Discrimination is the ability of a test to distinguish between different samples (genotypes)**

**Power of Discrimination**



**Speed of Analysis**

## Statistics of Chance: M&M Locus

### 6 Possible Alleles:

- Green
- Red
- Yellow
- Blue
- Brown
- Orange



## Probabilities



- **One allele from each parent means 2 copies of gene/locus**

$$\frac{1}{6} \times \frac{1}{6} = \frac{1}{36}$$

Frequency of any M&M genotype



## Probabilities

**Jolly Rancher:  
5 alleles**



**Mike & Ikes:  
5 alleles**




### Locus

M&M 6 alleles	$\frac{1}{6}$	$\times$	$\frac{1}{6}$	$=$	$\frac{1}{36}$
Jolly Rancher 5 alleles	$\frac{1}{5}$	$\times$	$\frac{1}{5}$	$=$	$\frac{1}{25}$
Mike & Ikes 5 alleles	$\frac{1}{5}$	$\times$	$\frac{1}{5}$	$=$	$\frac{1}{25}$

Chance an individual  
has a given genotype  $= \frac{1}{22,500}$

Chance 2 people have  
the same genotype  $= \frac{1}{506,250,000}$


**Who can't we  
exclude from  
the pool of  
suspects?**



[www.mms.com](http://www.mms.com)

© 2001 Mars, Incorporated. All Rights Reserved.

**The Power of  
Discrimination  
increases with  
the number of  
loci profiled**



$$0.044 \times 0.099 \times 0.050 = 2.18 \times 10^{-4}$$

$$1 / 2.18 \times 10^{-4}$$

**or 1 in 4591**

**The Power of Discrimination increases with the number of loci profiled**

<b>TPOX 8-12</b>	<b>0.044</b>
<b>D3S1358 16-17</b>	<b>0.099</b>
<b>FGA 21-23</b>	<b>0.050</b>
<hr/>	
<b>VWA 14-14</b>	<b>0.0088</b>
<hr/>	

**Random Match Probability**

**= 1 in  $5.3 \times 10^5$**

## Real-World Probabilities


- **Forensics labs use 13 different loci with multiple alleles**
- **Allele frequencies DO NOT follow mathematical principles - allele frequencies vary by population.**
- **These 13 loci allow for discrimination of any two people in the world (with the exception of identical twins), living or dead.**
- **Probability of a random match when all 13 loci typed: ~1 in 3 trillion.**

## TH01 Published Allele Frequencies by Population

<b>Allele</b>	<b>Caucasians n=302</b>	<b>African American n=258</b>	<b>Latinos n=140</b>
<b>5</b>	<b>.002</b>	<b>.004</b>	
<b>6</b>	<b>.232</b>	<b>.124</b>	<b>.214</b>
<b>7</b>	<b>.190</b>	<b>.421</b>	<b>.279</b>
<b>8</b>	<b>.084</b>	<b>.194</b>	<b>.096</b>
<b>9</b>	<b>.114</b>	<b>.151</b>	<b>.150</b>
<b>10</b>	<b>.008</b>	<b>.002</b>	<b>.014</b>
<b>11</b>	<b>.002</b>		

# **CODIS** **CO**mbined **DNA** **I**ndex **S**ystem

**A federally  
maintained  
database used  
by law  
enforcement  
officials**





**13 loci guarantees high power of  
discrimination**

## Real STR analysis

Four different fluorescent tags have been used to identify 7 amplified loci

Allele ladders are indicated by arrows





## Analysis of Results:

## Who can't be excluded?

