

Supplier Product Number

Shelf Life / Expiration / Use by Date

BDH83639.500 | Size: 5 L

HiPerSolv CHROMANORM

Lot: SAMPLE | Use by: 01.2015
Store at / Entreposage: Ambient Temp.
Date opened / Ouvert:

Acetonitrile

Acétonitrile
Acetonitrilo
for HPLC - SUPER GRADIENT Reag. Ph. Eur., ACS
water < 30 ppm - suitable for UPLC/UHPLC instruments

CAS: 75-05-8 | MW: 41.05 | Formula: C2H3N


Assay	Min. 99.9 %
Solution (100 g/l)	Neutral
Acidity	Max. 0.0008 meq/g
Alkalinity	Max. 0.0006 meq/g
Boiling point	90 → 82 °C
Colouration	Max. 10 APHA
Density (20°C)	0.781 → 0.784
n _D 20/D	1.34 → 1.345
Evaporation residue	Max. 2 %
Water	Max. 30 ppm
Absorbance (200 nm)	Max. 0.10
Transmittance (190 nm)	Min. 10 %
Transmittance (195 nm)	Min. 80 %
Transmittance (200 nm)	Min. 95 %
Transmittance (210 nm)	Min. 96 %

UN1648 ACETONITRILE



For Laboratory Use Only. Not for Drug, Food, or Household Use. | Pour usage en laboratoire uniquement. Non à usage pharmaceutique, alimentaire ou ménager. | Para uso en laboratorio únicamente. No para el uso de drogas, alimentos o casa.

For Specifications, C of A and SDS please go to vwr.com or scan the QR code on left. Pour spécifications, C of A et SDS s'il vous plaît aller à vwr.com ou scanner le code QR à gauche.



VWR ANALYTICAL



Danger - Highly flammable liquid and vapour. Harmful if swallowed, in contact with skin or if inhaled. Causes serious eye irritation. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician. Dispose of contents/container in accordance with local, state and federal regulations.

Danger - Liquide et vapeurs très inflammables. Nocif en cas d'ingestion, de contact cutané ou d'inhalation. Provoque une sévère irritation des yeux. Tenir à l'écart de la chaleur, des surfaces chaudes, des étincelles, des flammes nues et de toute autre source d'inflammation. Ne pas fumer. Porter des gants de protection/des vêtements de protection/un équipement de protection des yeux/du visage. EN CAS DE CONTACT AVEC LES YEUX: rincer avec précaution à l'eau pendant plusieurs minutes. Enlever les lentilles de contact si la victime en porte et si elles peuvent être facilement enlevées. Continuer à rincer. EN CAS d'exposition ou d'un malaise: Appeler immédiatement un CENTRE ANTIPOISON ou un médecin. Éliminer le contenu/récipient conformément aux règlements les, provinciaux et fédéraux.

Peligro - Líquido y vapores muy inflamables. Nocivo en caso de ingestión, contacto con la piel o inhalación. Provoca irritación ocular grave. Mantener alejado del calor, de superficies calientes, de chispas, de flamas abiertas y de cualquier otra fuente de ignición. No fumar. Llevar guantes/prendas/gafas/máscara de protección. EN CASO DE CONTACTO CON LOS OJOS: Aclarar cuidadosamente con agua durante varios minutos. Quitar las lentes de contacto, si lleva y resulta fácil. Seguir aclarando. EN CASO DE exposición o malestar: Llamar inmediatamente a un CENTRO DE INFORMACION TOXICOLÓGICA o a un médico. Deseche el contenido/envase conforme a las regulaciones locales, estatales y federales.

H+P numbers: +H225-H302+H312+H332-H319 -P210-P280-P305+P351+P338-P309+P310
Manufactured for: VWR International, LLC, 100 Matsonford Rd., Building One, Suite 200, Radnor, PA 19087, USA & VWR International Co, 2360 Argentinia Rd., Mississauga, ON, L5N 5Z7, Canada.
EMERGENCY/URGENCE: US (CHEMTREC): 800-424-9300; CA (CANUTEC): 613-996-6666
Made in France

0009

QR Code for Instant Access to Certificates and Safety Data Sheets

Lot/Batch Number